VERMEIREN

LE controller

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI

EN	Contents	
1.	Use	2
	Turn wheelchair On/Off	3
	Sleep mode	3
	Using the battery charge indicator.....	3
	Battery charging	4
	Adjusting the driving speed	5
	Using the horn	5
	The status indicator	5
	Using the joystick	5
2.	Troubleshooting	6

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
OPERATOR CONTROL LiNX LE

1. Use

The operator control built into your electric wheelchair enables you to control all the driving, steering, braking and operating processes. The wheelchair's electrical unit and electronics are constantly being monitored internally. Any fault in the electronics is indicated on the status indicator (1) (see the chapter on troubleshooting).

Vermeiren is responsible for the changes in the software. For changes in the software contact Vermeiren.

⚠ WARNING: Risk of burns – Be aware that the surface of the operator control can potentially get hot when it is exposed to strong sunlight for long periods.

Press the "ON/OFF" button (1). The charge indicator (2), which also indicates the batteries' charge level, lights up for a moment. If all the lights come on, the batteries are sufficiently charged. If the lights do not all come on, you must adjust your driving activities to suit the low capacity of the batteries. If the battery charger display lights do not come on, check the plug connections of the electronic system.

Now select the driving program you want by turning the speed switch (3), left turn to go slower and right turn to go faster.

To set the wheelchair in motion, push the joystick forwards slowly. The speed will increase when you push the joystick further forwards. If you want to turn left or right, simply push the joystick in the desired direction. To reverse, pull the joystick backwards from the central position.

⚠ Note that, when pressing the "ON/OFF" button, the joystick stays in the central position for at least five seconds. It is programmed this way for technical and safety reasons, to prevent the joystick being pushed at the same time as the wheelchair is being switched on. If both processes are activated at the same time, the steering unit locks itself and can only become operational again after being switched off completely.


You can sound an acoustic warning signal with the "Horn" (4) button.

To brake while driving, simply push the joystick to its central position. Move the joystick slowly for gradual braking. Simply release the lever for a quick stop; the wheelchair will then stop as quickly as possible.

Turn wheelchair On/Off

Turn the power on:


Press the On/Off button. The On/Off button will light green and all LEDs of the battery indicator turn on and display the current battery status.

If there is a fault with the system when powering up, the status indicator will indicate the fault with a series of red flashes. If the fault is one that prevents the system from driving, then the battery gauge will flash continuously.

If the wheelchair is turned on while the joystick is out of neutral, the wheelchair will not drive. Release the joystick back to neutral in 5 seconds and the fault will disappear. If the joystick is not released back to neutral (centre) within 5 seconds, the fault becomes a Module fault (the status indicator is flashing red), and the wheelchair will not drive. To clear this fault, loosen the joystick to the neutral position and turn the operator control off and then back on.

Turn the power off:

Press the On/Off button. The system will power down and the status indicator will switch off.

Emergency stop:

If the user needs to stop the wheelchair quickly, the On/Off button can be pressed to perform an EMERGENCY STOP. The wheelchair will come to a halt quickly.

Sleep mode


The operator control will go to sleep after a period of inactivity of a quarter.

When the system sleeps, it is partially turned off to reduce energy consumption.

To reactivate the operator control:

Press any button on the operator control or move the joystick. The system turns on.

Using the battery charge indicator

The battery charge indicator (located above the horn button) is used to indicate that your wheelchair power is on and provides an estimate of the remaining battery capacity.

Any green LED's indicate that your batteries are well charged.


If only amber and red LED's are light, the batteries are moderately charged. Recharge your batteries before you make a long trip.

If only red LED's are light, the batteries are running out of charge. Recharge as soon as possible.

The following table indicates what the battery charge indicator will display.

Display	Description	Meaning	Remarks
	All LEDs dark.	System off.	
	All LEDs are lit and the green LEDs are flashing.	System on. High voltage warning: battery voltage level has risen above the high voltage warning set-point.	Warning: battery too full. Slow down if you are driving down a slope.
	All LEDs are lit.	System on.	LEDs show the charge level.: Full battery (between 81 and 100% of full capacity)
	1 Green led is lit.	System on.	Battery level is decreasing.
	Only red and amber LEDs lit.	Low battery level.	Battery half full. Start the return journey.
	Red and amber LEDs lit.	Very low battery level.	Battery low (between 0 and 20% of full capacity). Recharge soon.
	Only Red LED lit.	Battery almost empty.	Recharge now.
	1 Red Led is flashing.	Battery empty. Low voltage warning: battery voltage level decreased below its low voltage warning set-point.	Recharge immediately.
	First red LED will flash on the battery indicator. The status indicator will display a flash code 2 or 7. The horn will sound once every 10 seconds.	Battery voltage decreases below the battery cut-off voltage.	Battery is fully discharged, and it is good practice to stop discharge at this point. The display and wheelchair shut down automatically.

Battery charging

Plug the battery charger into the remote's XLR socket.


The battery charger indicator will indicate the system is being charged by cycling between a left to right chase sequence, and then displaying the approximate battery charge state at the end of the chase sequence.

Driving is prevented (inhibited) while the system is being charged.

The LE system does not have to be powered up when charging the battery, however, if it is not powered up, then the battery gauge will not display the charging state / chase sequence.

Adjusting the driving speed

The speed switch allows the user to limit the maximum speed of the wheelchair (that is the speed when the joystick is fully deflected) their preference and environment.

The speed switch offers 10 discrete steps between the lowest speed (switch to the left) and the highest speed (switch to the right).


As a visual reminder, a speed symbol is positioned just below the speed switch to indicate the low and high positions.

Using the horn

Press the horn button (located above the ON/OFF button). The horn shall sound as long as the button is pressed.

The status indicator

The status indicator is located below the ON/OFF button. When the LiNX LE system is not powered up, the status indicator is not lit.

When the LiNX LE system is powered up, and there are no faults with the system, the status indicator will be lit green.

If, when powered up, there is a fault with the system, then the status indicator will flash red. The number of flashes will indicate the type of error. For flash codes refer to paragraph 2: Troubleshooting.


Using the Joystick

Moving the joystick will drive the electric wheelchair in that direction. The amount of joystick movement will determine the speed that the powerchair will move in that direction.

Make sure that the joystick is in the neutral (central) position when pressing the on/off button, otherwise the electronics will be locked. To take this block away let go the joystick until he goes back in neutral position.

Drive inhibit indication


Drive inhibit mode is indicated by the battery charger indicator with a right-to-left chase sequence.

The chase sequence starts with the green LED on the right-hand side, and by one-by-one, each LED will switch on and then off. When the sequence completes at the left-most red LED, it begins again at the righthand side.

The chase sequence continues until the error condition has been cleared.

2. Troubleshooting

The following diagnoses refer to all the electronics. The various modules for steering, power have their own light signals indicating the status of the module concerned (steady light = OK; flashing light = defective). The flashing code is read from the status indicator on the operator control.

If, when powered up, there is an error with the system, then the status indicator will flash red. The number of flashes will indicate the type of error. These is described in the table below.

This fault analysis aids problem analysis and problem solving by indicating the possible faults. If faults arise which cannot be identified and remedied with the aid of this analysis, kindly contact the manufacturer. We wish to point out that in any case of improper handling or if the delivered configuration has been changed, we cannot guarantee the accuracy of the fault indications.

Number	Problem/Fault	Checks
1	Remote / joystick error	1. Check that the joystick is placed in the neutral position when switch on the wheelchair. 2. Check that the joystick is not damaged and function very well.
2	Network or configuration error	1. Check all connection points, cables. 2. Verify software. 3. If you cannot found any fault contact your specialist dealer.
3	Left motor (or its connection) is defective.	1. Check whether the plugs of <u>both</u> motors have been inserted properly. 2. Check the plug contacts of <u>both</u> motors for corrosion or damage. 3. Check <u>both</u> motors. Unplug the motors and measure the plug connection with an Ohm meter. If you obtain readings of more than 1 Ohm or less than 100 milliOhms, the motor is defective. 4. Check the resistance of the motor. Measure every contact using an Ohm meter. ⚠ A Fault is also indicated if the motors are disconnected. Connect again and restart the system.
4	Right motor (or its connection) is defective.	As described above.
5	Left parking brake (or its connection) is defective or neutral active.	1. Check whether the motors' plugs have been inserted properly. 2. Check the plugs for corrosion or damage. 3. Check the parking brakes. Measure the resistance of the connections by means of an Ohm meter. If the resistance value is above 100 Ohms or beneath 20 Ohms, the parking brake is probably defective.
6	Right parking brake (or its connection) is defective.	As described above.
7	Module error (other than Remote)	One of the wheelchair parts can be damaged except the operator control itself.

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Vermeiren Suisse S.A.
Vermeiren Nederland B.V.