

eris



**VERMEIREN**  
we care for you

# eris

The Eris scooter is your in- and outdoor scooter for short distances to the shop, friends or family. Complete with indicators and lighting in LED version. The shopping tray can be taken off so you can carry it with you when purchasing small items at your local grocery store.



**DISPLAY**  
Easy to use display



**SHOPPING BASKET**  
Detachable



**STEERING COLUMN**  
Adjustable for greater driving comfort



**LED LIGHTING SYSTEM**  
Front and rear LED flashing light system with direction indicators



**SUSPENSION**  
Integrated shock absorption for a smooth driving experience



**PIVOTING SEAT**  
Equipped with a large pocket at the back

MAX. SPEED

Up to  
10 km/h

DRIVING RANGE

Up to  
37 km

USER WEIGHT

Up to  
140 kg



## TECHNICAL INFORMATION

Total width	570 mm
Total length	1185 mm
Total height	1150 mm
Seat height	450-500 mm
Turning radius	2900 mm
Driving range*	30-37 km
Motor	320 W
Speed*	6 or 10 km/h
Weight	75 kg
Max. user weight	140 kg
Batteries	2x12V - 38Ah/50Ah

\* The driving range is calculated according to European standards.

\* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1,5 kg and 1,5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



Contact your local supplier for more information.

## VERMEIREN GROUP

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



# VERMEIREN

we care for you