

Normandie XXL

Geriatric chair



VERMEIREN
we care for you

Optimal comfort

Care chair with adjustable backrest and legrest. Possible as a fixed chair or for patient transport with four optional casters. The height adjustable armrests can be fully lowered for lateral transfers while the anatomically shaped and foam padded seat and backrest allow for fatigue-free sitting and resting.

Options




Headrest



Extra wide armrests

Colours



Technical data

Total width	890 mm
Total height	1205 mm
Total length	820 mm
Seat width	730 mm
Seat height	510 mm
Seat angle	13°/40°

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Normandie XXL** complies with the highest norms, laid down by certified test institutes in most European countries.



Normandie XXL with wheels

Seat depth	490 mm
Height backrest	780 mm
Armrest height	0-220 mm
Weight	37,75 kg
Maximum user weight	350 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer