




Glue N° 11 into N° 10
with Anabond 202 instant
adhesive (3900649)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP

2991382

Inovys 2 Back Frame

Plan No.

Page 1/2

E000736/N

Scale: 1:5	Drafted	Verified	Confirmed
Date	22-06-2011	08-06-2026	08-06-2026
Initials	SD	MP	SD



Material			
No	Plan N°	Article N°	Qty
1	X001016	1001682	2
2	X001052	G066600	2
3	X001119	1001554	2
4	X003646	1999963	2
5	X004833	1745032	4
*	X004833	1745099	4
6	X004832	1744932	4
*	X004832	1744999	4
7	X004822	1744732	2
*	X004822	1744799	2
8	X004823	1744832	2
*	X004823	1744899	2
9	X004861	1906454	2
10	Y004353	F000411	2
11	Y017806	17562.. (40)	1
	Y017651	17560.. (45)	1
	Y017807	17563.. (50)	1
	YP06467	27347.. (35)	1
12	Y018198	1765462 (40)	1
	Y018199	1765562 (45)	1
	Y018200	1765662 (50)	1
	YP06468	2734862 (35)	1
13	M8x20	1001580	2
14	M8 DIN125	1003199	2
16	3,9x16 DIN7982	1005767	4
17	M6x25 ISO7380	1006998	2
18	M6x25 DIN912	1001494	4
*	M6x25 DIN912 Black	2996887	4
19	M6x20	1009932	2
20	Ø22	1990648	4
21	D8x40	1992510	1
22	X005513	1992600	2
24	---	1996806	2

(40 - 45 - 50 - 35 = Seat width)

Material			
No	Plan N°	Article N°	Qty
25	M6x16 DIN7991	1994752	4
*	M6x16 DIN7991 Black	2994064	4
28	M8x20 ISO7380	1992349	2
29	L34	1906938 (40)	1
	L34	1906940 (45)	1
	L34	1906942 (50)	1
	L34	2996478 (35)	1

* - Use for Inovys Junior (35)

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP

Inovys 2 Back Frame Page 2/2				2991382
				Plan No.
				E000736/N
Scale: 1:5	Drafted	Verified	Confirmed	
Date	22-06-2011	08-06-2026	08-06-2026	
Initials	SD	MP	SD	