




INOVYS 2 LOWER FRAME 24 INCH

Scale: 1:6	Drafted	Verified	Confirmed
Date	27.09.2011	08/06/2026	08/06/2026
Initials	SD	MP	SD

Sheet 1/2

Plan No.

E000747/M



Material			
No	Plan N°	Article N°	Qty
1	Y010736	16386.. (40)	1
	Y011927	16596.. (45)	1
	Y011934	16597.. (50)	1
	YP06450	27342.. (35)	1
2	Y004353	F000411	2
3	Y011640	1655154	1
4	Y011652	1655254	1
5	Y017901	1758901	2
6	X000188	G078500	2
7	X001380	1009726	8 (a) 4 (b)
8	X001016	1001682	1
9	---	1905244	2
10	X002357	1993126	2 (b)
11	X002437	1992699	2
12	X003451	1658800	2
13	X003452	1658900	2
14	X003596	1998528	2
15	X005720	1905243	2
16	---	1000181	2
17	M12 DIN125	1001561	2
18	---	1001610	2
19	M6 DIN984	1001680	1
20	M12 DIN984	1001684	2
21	M6X35 ISO7380	1006057	1
22	M12X75 DIN931	1009991	2
23	M7	1900008	1
24	200x50	1914609	2
25	6001-2Z	1992559	4
26	---	1992599	2
27	M16 type S DIN17222	1911653	2
28	M8x45 ISO7380	1994473	1
29	M8X76 DIN912	1998993	2
30	M8	1998994	2

Material			
No	Plan N°	Article N°	Qty
31	24"	1914602	2
32	400N	1005602	1
33	Left	1900915	1
34	Right	1900916	1
35	X003800	1900347	2

(a) For wheel in rear position
(b) For wheel in front position

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.				Sheet 2/2	---
INOVYS 2 LOWER FRAME 24 INCH				Plan No.	E000747/M
Scale: 1:6	Drafted	Verified	Confirmed		
Date	27.09.2011	08/06/2026	08/06/2026		
Initials	SD	MP	SD		