



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

2993366

Plan No.

E001056/I

MOTOR SET INOVYS 2-E

Scale:	11:50	Drafted	Verified	Confirmed
Date	14/10/2016	24/10/2025	24/10/2025	
Initials	SD	PM	AS	



Material			
No.	Plan	Article	Quantity
1	Y021572	2022954	1
2	Y021954	2023154	1
3	Y021955	2022854	2
4	M4 DIN985	2994901	2
5	Y022438	2049354	1
6	M5 DIN985	2994006	6
7	M6 DIN985	2994063	4
8	M8 DIN985	1996213	4
9**	X005710	2026200	2
10	Y022608	F007420	2
11	M4x12 DIN963	1919449	2
12	M6x16 ISO7380	2994067	6
13	∅8x40	1992510	1
14	M6x30 ISO7380	2994129	2
15	M5x30 DIN912	1919450	4
16	Act. TA9	1912783	2
17	Batt. TBB2	1912784	1
18*	Charg. TP8-11A3	1914468	1
19	Hand. TH21	1912786	1
20	Hub TJB5	1912787	1
21	Cable TYC	1914469	1
22	M8x50 ISO7380	1915952	1
23	M5x20 DIN912	1915452	2
24	Y022168	2038754	1
25	∅4x16	1919451	2
26*	---	1912887	1
27	X005575	2023600	2
28	M8x55 ISO7380	1919452	1
29	M8x45 ISO7380	1915945	1
30**	M8x80 ISO7380	1919453	2
31**	Y022603	2026354	1
32**	Y022604	2026454	1
33**	M12	1916509	2
34**	M8 DIN126	1996214	2

* Not depicted

** Used for seat width 40 cm, for SW 45 and 50 see E000383

Item 35: see E000383 for article numbers for SW 40-45-50

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

2993366

MOTOR SET INOVYS 2-E

Plan No.

E001056/I

Scale: 11:50	Drafted	Verified	Confirmed
Date	14/10/2016	24/10/2025	24/10/2025
Initials	SD	PM	AS

