



Inovys II

L70

Granulated seating system

VERMEIREN created the L70 positioning system in combination with the existing Inovys II wheelchair. For people who are bound to a passive sitting position, the Inovys L70 offers ideal seating comfort.

The polystyrene-filled cushions in the back, seat, and leg area ensure pressure distribution and follow the contours of the body.



VERMEIREN
we care for you

Pleasant seating comfort for passive wheelchair users.

VERMEIREN can offer you a complete and highly dynamic solution that provides maximum support for the entire body. All cushion elements are provided with a zipper so that they can be easily removed and cleaned.

The **Inovys II L70** is also equipped with a fully adjustable headrest both in height and in angle and this without additional tools. Optionally, an L58G can also be mounted that automatically assumes any desired shape. A range of options means that the Inovys II L70 can be fully adapted to the needs and wishes of almost every user.

Individual adjustments are perfectly possible due to the adjustability of the angle of the foot section, the leg section, the seat section and the inclination of the back section. The drum brakes for the attendant and additional parking brakes ensure that the user can sit in the chair in complete safety.

An optional **washable cover** in gray terry cloth is available and thus ensures more hygiene during daily use of the Inovys II L70.



⊕ All cushions are filled with resilient polystyrene granules that provide perfect support for the body and connect to the its contours.

⊕ By adding or reducing the polystyrene granules, the pressure on body contours can be increased or decreased.

⊕ Due to the various adjustment options (head, back, seat, and leg section), the wheelchair can be optimally adjusted to each user.

⊕ The comfort arms provide a secure and pleasant feeling that allows the user to completely relax.



Change to a personalised resting position thanks to the independent adjustment possibilities of seat, back rest, feet and legs.

3 independent granulated compartments



Each compartment features bead-filled cushions that provide optimum body & contour support.
- Chest compartment -



- Seat compartment -



- Leg compartment -



Easy removal of the compartments by means of a zipper



B12 - Table



B20 - Safety belt



B52 - Serum holder



B54 - Cloth



B74 - Drum brakes



B77 - Anti-tipping



B83 - Canister holder



L58G - Comfort headrest



L71 - Granulated removable armrests

options

Granulated filling



To offer a maximum seating comfort **VERMEIREN** has chosen to select a premium quality of granules that reach a density up to 15 kg per square meter. Each granule is between 3,6 mm and 5,6 mm in diameter and has an insulating effect creating a comfortable feeling when in contact with a human body. These granules are manufactured from 100% expanded polystyrene and have been certified by the National Institute of Hygiene that they do not pose any hazard to human health. These granules can be recycled and as such, leave a minimal ecological footprint.



Refill : After a longer usage you may feel that your pillows need additional filling or even lesser depending on your needs. Contact your trusted **VERMEIREN** dealer and ask for granules to refill your cushions and retighten your positioning.

Full circle

If it needs to move, it needs wheels. And wheels are one of the key elements of our products. One solid axle combined with spokes, rims and a tube filled with air. That's all it takes to move us forward. Very basic but ever so essential. And that is the reason why **VERMEIREN** decided to produce this key element internally.

From small front wheels to big rear wheels to solid rubber wheels for power wheelchairs... Every part is manufactured internally so that we can keep a close eye on the production and even more on the consistent quality that do us proud.

Technical data

Total width	680 - 730 mm
Total height (without headrest)	990 - 1060 mm
Total length	1280 mm
Seat inclination	0° - 30°
Legrest - footplate inclination	73° - 107°
Height armrests	220 - 320 mm



Insulated material provides a warm and comfortable feeling when seated.


Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Inovys II L70** complies with the highest norms, laid down by certified test institutes in most European countries.

Seat height	510 - 580 mm
Height backrest	550 - 600 mm
Seat depth	460 - 520 mm
Weight	38 kg
Maximum user weight	135 kg
Additional info	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you