



GEMINI 2

Quality and possibilities.

The beauty of comfort and simplicity.

The completely adjustable and customizable
Gemini 2 is a perfect companion for family
outings or in-home interaction.



Easy & Smart positioning control

■ A correct posture, even at a very young age, is important.

Thanks to our wide range of positioning accessories the **Gemini 2** can be modified so that the experience becomes personal.

Possible adjustments include:

- Head corrections
- Trunk corrections
- Hip corrections
- Feet corrections



Efficient postural change

- For short or long usage, postural change through the day is a must for your kid. Either for pressure relief, blood circulation or just for comfort.

Whether you like to be in permanent visual contact with your kid, or your kid prefers to see and explore the world in front of him: easily reverse the seat with the smart **ReverSeat** system.



Make your own Gemini

- Completely personalize and adjust the **Gemini 2** with a wide range of options and accessories.



Table - U90



Rain cover - U96

Upholsteries & Frame colours



Black & Grey



Blue & Grey



Green & Grey



Purple & Grey



Blue & Black 3D Net



Frame colour: Matt Black



Frame colour: Dotted Dark Grey

Easy to use

- Configure and easily adapt your **Gemini 2** with the smart tool-free system. Pull, twist or turn the knobs and reconfigure your chair in seconds.

The **Gemini 2** has been redesigned with compactness in mind, providing easy storage and transport.



Pull and reverse



Backrest adjustment



Pull and fold



Footpedal brake



Adjustable pushbar



Adjustable footplate



Ingenious 2-in-1 system

- Create 2 chairs with 1 seat system: an indoor chair and a stroller. Using the hi-low indoor frame, you can quickly mount the seat of Gemini 2 on top. Turn your stroller into an inclinable chair for many fun activities at the dining table.





8" Inflatable/PU front wheels



12" Inflatable/PU rear wheels



Hip cushion 2 cm - B28A
Hip cushion 4 cm - B28B



Trunk cushion 2 cm - B28C
Trunk cushion 4 cm - B28D



Abduction block - U53



Shopping basket - U61



Legs cover - U62



Canopy - U67



Directional wheel lock - U68



Barrier - U64

180° adjustable
pushing bar

Foldable and
removable canopy - U67



Adjustable
head support - U72

Double-sided
safety harness


Side Support
(hips/trunk) - U69

Tubular armrests

PU or inflatable rear wheels

SCAN ME to EXPLORE ME





Respirator shell - U70



Trunk stabilization profile - U71



Tubular arm rest - U84




Adjustable trunk pelots - U85



High-low indoor frame - U82



Adjustable hip pelots - U86



Additional small trunk pelots - U87



Frame protection - U88



Lumbar cushion - U89



Plastic table - U90



Feet positioning system - U91



Technical data Gemini 2 (for models with seat width 320 and 400 mm)

| | |
|-------------------|-------------------|
| ■ Total width | 620 mm & 690 mm |
| ■ Total height | 1200 mm & 1275 mm |
| ■ Total length | 950 mm & 1060 mm |
| ■ Folded width | 620 mm & 690 mm |
| ■ Seat width | 320 mm & 400 mm |
| ■ Height armrests | 200 mm |

| | |
|-----------------------|-----------------------------|
| ■ Seat height | 510 mm |
| ■ Height backrest | 520 - 640 mm & 650 - 770 mm |
| ■ Seat depth | 300 - 360 mm & 360 - 420 mm |
| ■ Weight | 18 kg & 27 kg |
| ■ Maximum user weight | 35 kg & 50 kg |
| ■ Additional info | CE |

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Gemini 2** complies with the highest norms, laid down by certified test institutes in most European countries.



VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer