

Squod s.u.

ELECTRONIC STANDUP WHEELCHAIR

The Squod S.U. is an electronic stand up wheelchair complete with lights both front and back, electric adjustable backrest (SE41) and an electronic stand up system with anti-tipping wheels and a footplate out of one piece in front. Comes standard with taxi hooks and two motors of 200W and 2 batteries of 75 Ah.

Compatible iPortal*



Adjustable back 30°



C55
inox



C30
diamond
grey



C41
green

model SQUOD



Upholstered seat and back



Shark 2
Controller

In- and outwards
foldable legrests

Foldable footplates

Comfort

Ideal compromise for indoor outdoor usage. The standing function will adapt to each situation.

The Squod is equipped with standard lighting.



Anti-tipping wheels



Armrests adjustable in
height and width



2x12V
75Ah



2x200W



6 km/h



30 km



6°/10,5%



60



1400



Power
Buddy

Wide range of accessories and options.

*iPortal is compatible with all the electronic wheelchairs of Vermeiren equipped with a DX2 Basic, DX2 LCD or special controller.



Squod

600
650
710

390-430
440-490
500-550

200-250

600

440

510

1070

111

130

Squod SU

650

450

min. 220
&
max. 300

580

min. 400
&
max. 450

510 or 610

1240

121

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

