

Timix s.u.

Compatible iPortal*



Fully collapsible knee-lock

ELECTRONIC STANDUP WHEELCHAIR

The Timix S.U., also known as the Timix Stand Up, is an electronic stand-up wheelchair completely controllable from its controller. From a vertical stand-up position to a complete horizontal position in just one movement.



C71
blue



C72
red



C73
green



Chest safety belt with velcro

Back and seat in soft version

Automatic length compensating footplate

Anti-tipping wheels



2x12V-80Ah	2x220W	6 km/h	63 km	9°	100	1320
	2x350W	10 km/h	40 km			
		14 km/h	39 km			



Wide range of accessories and options.

*iPortal is compatible with all the electronic wheelchairs of Vermeiren equipped with a DX2 Basic, DX2 LCD or special controller.



670

450

min. 242
&
max. 350

540
(without cushion)

450
500

610

1260

160

130

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.



TIMIX S.U.