

Eclips X4

Multi adjustable manual wheelchair



Eclips X4 is a VERMEIREN classic wheelchair that has proven its quality throughout history. It is a light aluminum wheelchair with a multi-adjustable frame. The seat is adjustable in depth, has a reinforced cross for extra strength and Quick-release rear wheels that will provide a speedy disassembly when placing the wheelchair into a car. The footrests can be folded both inward and outward and can be removed to limit the weight of the heaviest part.




VERMEIREN
we care for you

Standard

- Adjustment of the seat depth
- Aluminum, height adjustable forks to keep the wheelchairs weight minimal
- Foldable footplates
- Rear wheels with Quick-release system
- Backrest adjustable in height

Colours



Features




Extra long and padded armrests



Depth adjustable frame



Height adjustable frontwheels and aluminium forks



Tipping piece at the rear to overcome obstacles

Technical data

■ Total width	590/620/640/660/680/700/740 mm
■ Total height	910 - 1015 mm
■ Total length	1000 - 1210 mm
■ Folded width	300 mm
■ Seat width	390/420/440/460/480/500/550 mm
■ Height armrests	220 mm

■ Seat height	390 - 520 mm
■ Height backrest	420 - 460 mm
■ Seat depth	420 - 500 mm
■ Weight	17 kg
■ Maximum user weight	150 kg
■ Additional info	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Eclips X4** complies with the highest norms, laid down by certified test institutes in most European countries.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you