

# one<sup>+</sup>

Compact indoor scooter  
Perfect for daily trips  
Adapted to your style



one<sup>+</sup> = 4

4 different colours included. Change your scooter's style whenever you want. No tool needed.




**VERMEIREN**  
we care for you

## Compact and stylish

**One+** is a perfect compact and stylish scooter for in- and outdoor use. Adapted to your style and your lifestyle! Lightweight and compact in its form, **One+** can be used for short daily trips to restaurants, stores or in small spaces. The scooter fully disassembles for compact travelling or storage. 4 unique colours can be interchanged easily for a matching scooter to your style, at any time.

## Colours

Your **One+** scooter comes with 4 interchangeable colours included. Easy to change and adjustable to style. One = 4



blue

green

orange

pink

## Products designed to suit you

**VERMEIREN** always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **One+** complies with the highest norms, laid down by certified test institutes in most European countries.



Fully dismantlable in a matter of seconds



User-friendly controller with a battery status LED indicator



Adjustable steering column



Fits in a small car

## Technical data

Total width	490 mm
Total height	880 mm
Total length	1030 mm
Seat height	525-570 mm
Turning radius	1150 mm
Battery	24V - 25 Ah

Driving range	30 km*
Motors	270 W
Speed	8 km/h*
Weight	50 kg
Maximum user weight	136 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

\* Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%).

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you