

The logo features a stylized blue 'C' with a white eye-like shape inside, followed by the text 'ECLIPS X2 KIDS' in a bold, blue, sans-serif font.

ECLIPS X2 KIDS

Lightweight and foldable wheelchair,
made in aluminium.





Compact and versatile

Compact and versatile, the **Eclips X2 Kids** comes packed with all elementary features to tweak this chair to the personal needs of its user.

Eclips X2 Kids is completely foldable for easy transport, and all of this comes in a trendy, eye-catching red colour.

Equipped with brakes, foldable backrest and safety belt. Available in seat width 32 and 36 cm.





Foldable backrest



Very compact



BZ7 - Optional legrest





Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Eclips X2 Kids** complies with the highest norms, laid down by certified test institutes in most European countries.

Technical data

■ Total width	520 mm
■ Total height	930 mm
■ Total length	830 mm
■ Folded width	320 mm
■ Seat width	320 - 360 mm
■ Height armrests	190 mm

■ Seat height	440 mm
■ Height backrest	330 mm
■ Seat depth	320 mm
■ Weight	13,6 kg
■ Maximum user weight	60 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you