

# Alesia

## Multifunctional geriatric chair

The multifunctional Alesia chair offers a reclinable backrest with a dual control. Adjustable seat inclination, adjustable armrests, retractable footplate and footbrakes at the rear. Available in a wide range of different colours and a variety of options to upgrade the chair.



**VERMEIREN**  
we care for you

## Alesia, total comfort

Alesia is fitted with 4 swivel wheels with a diameter of 125 mm to facilitate moving around.

Alesia exists in different versions: with manually adjustable legrests or without legrests. The footplate can be retracted inwards and prevents tipping over. The armrests are adjustable in height and completely retractable facilitating transfers.

Complete with centralised brakes, rear directional wheels and a reinforced back which can be inclined by the user or by a third person (the controllers are interconnected).

## Colours



## Options

- AL04 - side-supports
- AL12 - table
- AL52 - serum holder
- AL55 - adjustable and removable headrest
- AL83 - canister holder
- AL84 - document holder
- AL85 - urine bag holder



## Technical data

Total width	700 mm
Total length	800 mm
Total height	1330 mm
Seat width	540 mm
Seat height	530 mm
Seat angle	12°

## Products designed to suit you

**VERMEIREN** always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, **Alesia** complies with the highest norms, laid down by certified test institutes in most European countries.



Seat depth	480 mm
Height backrest	800 mm
Height armrests	0 - 270 mm
Weight	33 kg
Maximum user weight	140 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer