

Navix (FWD)

FEATURES

Better manoeuvring with a compact front-wheel drive

ELECTRONIC ADULT WHEELCHAIR

The Navix allows for sharper turning indoor and for better overcoming obstacles like thresholds when entering the local grocery store. The user friendly Navix is adjustable to individual needs and allows for easy transfers to bed, bath, toilet, car etc..

C30 metallic grey	C71 blue	C72 red	C73 green

REARWHEEL DRIVE available



- Better overcoming obstacles
- Better taking sharp turns
- Easy to reassemble
- Total width for Navix limited to 570 mm
- Ground clearance of 60 mm
- Adjustable to your needs



Optional electrical inclination of the seat and backrest



Optional 300 mm lift for a better accessibility



Optional: BZ7-E electrical adjustable legrests

Better grips on outdoor surfaces with profiled large front wheels

2x12V 38Ah	2x220W	6 km/h	20 km	9°	50	1250



--	--	--	--	--	--	--	--	--	--

580	400	260-350	400/450/500	380-500	580	330-400	1010 (945 RWA)	79 kg	130
580	450	260-350	400/450/500	380-500	580	330-400	1010 (945 RWA)	79 kg	130
630	500	260-350	400/450/500	380-500	580	330-400	1010 (945 RWA)	79 kg	130



RWD: Rear wheel drive - Sizes in mm, weight in kg. Deviations in sizes and weight (1,5 mm, 1,5 kg and 1,5°) are possible.