



Material			
No.	Plan	Article	Quantity
1	X001016	1001682	2
2	X001119	1001554	4
3	Y018754	17769..	1
4	Ø16	1001598	4
5	Ø19	1001602	2
6	M8x25 DIN912/8.8	1001499	2
7	M8 DIN127	1001570	2
8	M8x12 DIN912	1008963	2



E000846 / D

LUGGAGE RACK TYPE TIMIX BAGAGE REK TYPE TIMIX				Plan No.
				E000846 / D
Scale: 1:2	Drafted	Verified	Confirmed	
Date	27/11/2012	19/12/2014	19/12/2014	
Initials	SD	TT	PJ	