



Material				
No.	Plan	Article	Quantity	
1	P00067819	2679604 (38)	1	
	P00067831	2680104 (39)	1	
	YP06392	2729504 (40)	1	
	P00067821	2679704 (42)	1	
	P00067823	2679804 (44)	1	
	P00067825	2679904 (46)	1	
	P00067827	2680004 (48)	1	
	YP06394	2729704 (50)	1	
	2	XP02084	GP16000	2
	3	X001052	G066600	2
4	X001473	1002994	2	
5	X004861	1906454	2	
6	M4x12 DIN963	1005401	4	
7	STAR KNOB M6x20	1009932	2	
8	M5 DIN985	2994006	2	
9	M6 DIN125-1	2994068	2	
10	M5x25 DIN7991	2994071	2	
11*	Velcro	1009792	1[m]	

*- Not depicted.

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

ALU BACK FRAME D200/V300				---
				Plan No.
				E00067829/A1.1
Scale:	1:3	Drafted	Verified	Confirmed
Date	2023-09-26	26/09/2023	26/09/2023	
Initials	AS	PM	AS	

