



DETAIL 11
SCALE 2X

* Item 21 add torque 60Nm

EC00066/M

SIDE FRAME TYPE JAZZ S50				Plan No. / Plannr.
				EC00066/M
1/2				
Scale / Schaal: 1 : 8	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	03/11/2017	24/02/2022	24/02/2022	
Init. / Init.	SL	AZ	JJQ	

Material			
NO.	Plan	Article	Number
1	YC00358(L)	1C158	1
	YC00359(R)	1C159	1
2	---	1908236	2
3	---	2990070	4
4	X002679	1994452	2
5	X000582	1005740	2
6	---	1008492	8
7	---	1001564	8
8	---	1006850	4
9	X004309	1904133	2
10	---	1991017	2
11	(L)	3990524	1
	(R)	3990525	1
11A	XC00119(R)	3990382	1
	XC00120(L)	3990383	1
11B	---	1992597	2
11C	X002163	1488100	2
11D	X002169	1992791	2
11E	X002132	1992423	2
16	X004415	17193	2
17	---	1991447	2
18	---	1001572	2

Material			
NO.	Plan	Article	Number
19	---	1998994	2
20	---	1998993	2
21	---	3994158	2
22			
23	X003293	1001626	2
24a	XC00122	3990395*	2
24b	XC00037	3990106*	2
25	---	1001684*	2
26	X003801	1900365*	2
27			
28	X004308	1904129	4
29	---	1903903	2

*If assembly rear wheel fixed axle, take out 24b;
If assembly rear wheel quick released axle,
take out 24a&25&26, only use 24b.

EC00066/M

SIDE FRAME TYPE JAZZ S50				Plan No. / Plannr.
				EC00066/M
2/2				
Scale / Schaal: 1 : 8	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	03/11/2017	24/02/2022	24/02/2022	
Init. / Init.	SL	AZ	JJQ	