



EP00614 / B

Isometric view
Scale: 1:10

BACKREST GREY/GREEN GEMINI II 32 2020: 2996830			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,75m
2	FABRIC ALCALA (63) GREEN	2996798	0,2m
3	NYLON BLACK	1008943	0,3m
4	VELCRO BLACK 25 MM	1009709	0,3m
5	VELCRO HOOK BLACK 25 MM	1009710	0,3m
6	XP03151	2996723	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,01rl
8	VULLING CARDAFILL	1991745	0,01kg
9	VELCRO BLACK 100MM	2991789	0,3m
10	VELCRO HOOK BLACK 50 MM	1009791	0,8m
11	VELCRO BLACK 50 MM	1009792	0,6m

BACKREST GREY/BLUE GEMINI II 32 2020 : 2996831			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,75m
2	FABRIC ALCALA (65) BLUE	2996797	0,2m
3	NYLON BLACK	1008943	0,3m
4	VELCRO BLACK 25 MM	1009709	0,3m
5	VELCRO HOOK BLACK 25 MM	1009710	0,3m
6	XP03151	2996723	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,01rl
8	VULLING CARDAFILL	1991745	0,01kg
9	VELCRO BLACK 100MM	2991789	0,3m
10	VELCRO HOOK BLACK 50 MM	1009791	0,8m
11	VELCRO BLACK 50 MM	1009792	0,6m

BACKREST GEMINI 2020				Plan No.
Page 1/3				EP00614 / B
Scale: 1:10	Drafted	Verified	Confirmed	
Date	21/09/2021	26/10/2021	26/10/2021	
Initials	OY	AZ	TK	

BACKREST GREY/BLACK GEMINI II 32 2020: 2996832			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,75m
2	FABRIC ALCALA (14) BLACK	2994281	0,2m
3	NYLON BLACK	1008943	0,3m
4	VELCRO BLACK 25 MM	1009709	0,3m
5	VELCRO HOOK BLACK 25 MM	1009710	0,3m
6	XP03151	2996723	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,01rl
8	VULLING CARDAFILL	1991745	0,01kg
9	VELCRO BLACK 100MM	2991789	0,3m
10	VELCRO HOOK BLACK 50 MM	1009791	0,8m
11	VELCRO BLACK 50 MM	1009792	0,6m

BACKREST BLACK/BLUE AIRFL GEMINI II 32 2020 : 2996834			
No.	Plan	Article	Quantity
1	FABRIC AIR FLOW BLACK 900	2992087	0,75m
2	FABRIC AIR FLOW BLUE 9570	2992088	0,2m
3	NYLON BLACK	1008943	0,3m
4	VELCRO BLACK 25 MM	1009709	0,3m
5	VELCRO HOOK BLACK 25 MM	1009710	0,3m
6	XP03151	2996723	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,01rl
8	VULLING CARDAFILL	1991745	0,01kg
9	VELCRO BLACK 100MM	2991789	0,3m
10	VELCRO HOOK BLACK 50 MM	1009791	0,8m
11	VELCRO BLACK 50 MM	1009792	0,6m

BACKREST GREY/PURPLE GEMINI II 32 2020 : 2996833			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,75m
2	FABRIC ALCALA (62) PURPLE	2996800	0,2m
3	NYLON BLACK	1008943	0,3m
4	VELCRO BLACK 25 MM	1009709	0,3m
5	VELCRO HOOK BLACK 25 MM	1009710	0,3m
6	XP03151	2996723	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,01rl
8	VULLING CARDAFILL	1991745	0,01kg
9	VELCRO BLACK 100MM	2991789	0,3m
10	VELCRO HOOK BLACK 50 MM	1009791	0,8m
11	VELCRO BLACK 50 MM	1009792	0,6m

BACKREST GREY/BLUE GEMINI II 40 2020 : 2996973			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,85m
2	FABRIC ALCALA (65) BLUE	2996797	0,25m
3	NYLON BLACK	1008943	0,4m
4	VELCRO BLACK 25 MM	1009709	0,4m
5	VELCRO HOOK BLACK 25 MM	1009710	0,4m
6	XP03179	2996796	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,02rl
8	VULLING CARDAFILL	1991745	0,02kg
9	VELCRO BLACK 100MM	2991789	0,4m
10	VELCRO HOOK BLACK 50 MM	1009791	1m
11	VELCRO BLACK 50 MM	1009792	0,8m

EP00614 / B

BACKREST GEMINI 2020				Plan No.
Page 2/3				EP00614 / B
Scale: 1:10	Drafted	Verified	Confirmed	
Date	21/09/2021	26/10/2021	26/10/2021	
Initials	OY	AZ	TK	

BACKREST GREY/BLACK GEMINI II 40 2020: 2996974			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,85m
2	FABRIC ALCALA (14) BLACK	2994281	0,25m
3	NYLON BLACK	1008943	0,4m
4	VELCRO BLACK 25 MM	1009709	0,4m
5	VELCRO HOOK BLACK 25 MM	1009710	0,4m
16	XP03179	2996796	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,02rl
8	VULLING CARDAFILL	1991745	0,02kg
9	VELCRO BLACK 100MM	2991789	0,4m
10	VELCRO HOOK BLACK 50 MM	1009791	1m
11	VELCRO BLACK 50 MM	1009792	0,8m

BACKREST BLACK/BLUE AIRFL GEMINI II 40 2020 : 2996976			
No.	Plan	Article	Quantity
1	FABRIC AIR FLOW BLACK 900	2992087	0,85m
2	FABRIC AIR FLOW BLUE 9570	2992088	0,25m
3	NYLON BLACK	1008943	0,4m
4	VELCRO BLACK 25 MM	1009709	0,4m
5	VELCRO HOOK BLACK 25 MM	1009710	0,4m
16	XP03179	2996796	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,02rl
8	VULLING CARDAFILL	1991745	0,02kg
9	VELCRO BLACK 100MM	2991789	0,4m
10	VELCRO HOOK BLACK 50 MM	1009791	1m
11	VELCRO BLACK 50 MM	1009792	0,8m

BACKREST GREY/PURPLE GEMINI II 40 2020 : 2996975			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,85m
2	FABRIC ALCALA (62) PURPLE	2996800	0,25m
3	NYLON BLACK	1008943	0,4m
4	VELCRO BLACK 25 MM	1009709	0,4m
5	VELCRO HOOK BLACK 25 MM	1009710	0,4m
16	XP03179	2996796	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,02rl
8	VULLING CARDAFILL	1991745	0,02kg
9	VELCRO BLACK 100MM	2991789	0,4m
10	VELCRO HOOK BLACK 50 MM	1009791	1m
11	VELCRO BLACK 50 MM	1009792	0,8m

BACKREST GREY/GREEN GEMINI II 40 2020: 2996972			
No.	Plan	Article	Quantity
1	FABRIC ALCALA (22) LIGHT GREY	2996799	0,85m
2	FABRIC ALCALA (63) GREEN	2996798	0,25m
3	NYLON BLACK	1008943	0,4m
4	VELCRO BLACK 25 MM	1009709	0,4m
5	VELCRO HOOK BLACK 25 MM	1009710	0,4m
16	XP03179	2996796	1
7	MOUSSE ON ROLL 5MM 180MX2,10M	1990296	0,02rl
8	VULLING CARDAFILL	1991745	0,02kg
9	VELCRO BLACK 100MM	2991789	0,4m
10	VELCRO HOOK BLACK 50 MM	1009791	1m
11	VELCRO BLACK 50 MM	1009792	0,8m

EP00614 / B

BACKREST GEMINI 2020				Plan No.
Page 3/3				EP00614 / B
Scale: 1:10	Drafted	Verified	Confirmed	
Date	21/09/2021	26/10/2021	26/10/2021	
Initials	OY	AZ	TK	