



LEGREST NORMANDIE XXL				Plan Number
				Sheet 1/2 E000607/B
Scale: 1:5	Drafted	Verified	Confirmed	
Date	03/12/2009	23/09/2022	23/09/2022	
Initials	DDR	PM	AZ	



Material			
No.	Plan	Article	Quantity
1	Legrest Blue	2997482	1
	Legrest Lavender	2997483	1
	Legrest Brown	2997484	1
	Legrest Black	2997485	1
	Legrest Grey	2997486	1
2*	Cover Blue	2997487	1
	Cover Lavender	2997488	1
	Cover Brown	2997489	1
	Cover Black	2997490	1
	Cover Grey	2997491	1
3*	X004507	1905152	1
4*	Insert M6x10	1992516	4
5*	X004502	H7621	1
6*	---	1005328	1
7	DIN125A D6,4	1001564	4
8	M6x16	1008459	4

*Subparts of complete legrest



LEGREST NORMANDIE XXL				Plan Number
				Sheet 2/2 E000607/B
Scale: 1:4	Drafted	Verified	Confirmed	
Date	03/12/2009	23/09/2022	23/09/2022	
Initials	DDR	PM	AZ	