



Drawing shows only right side

EP00378 / B

ADJUSTABLE BACKREST WITH FOLDABLE HANDLE TRIGO				Plan No. EP00378/B
Page1/2				
Scale: 1:5	Drafted	Verified	Confirmed	
Date	14/02/2019	10/12/2019	10/12/2019	
Initials	AS	FD	TK	

Material			
No.	Plan	Article	Quantity
1	YP06265 (25-35)	2785265	1
	YP06268 (35-45)	2785365	1
	YP06269 (45-50)	2785465	1
** 2	YP07062(32-37)	2744565	1
	YP07060(35-40)	2744665	1
	YP07061(40-45)	2744765	1
	YP07063(45-50)	2744865	1
3	M6x30	1992721	1
4	Washer M6	2994068	1
5	M6	2994059	3
6	XP01560	2993561	4
7	M6x40	2994066	2
8	Fold. push handle	EP00246	1
9	Y004353	F000411	1
10	XP02043	2993671	1
11	M4	1001678	1
12	M4x16	2994070	1
13	M6x10	2994127	3
14	M6X10	2994284	2
* 15	Cap D22	1990648	1

\* When unused folding handle / Gdy nie uzywana raczka skladana.

\*\* Gdy nie stosujemy poz. 8 to nalezy uzyc detalu o rozmiar wyzej.

\*\*If no. 8 isn't used, take detail one size bigger.

EP00378 / B

ADJUSTABLE BACKREST WITH FOLDABLE HANDLE TRIGO				Plan No.
				Page 2/2 EP00378/ B
Scale: 1:5	Drafted	Verified	Confirmed	
Date	14/02/2019	10/12/2019	10/12/2019	
Initials	AS	FD	TK	