



E000566 / F

MOTOR 6KM/h SQUOD / SQOUD 2009 / SQOUD SU				Plan No.
				E000566 / F
Scale: 1:5	Drafted	Verified	Confirmed	
Date	11/02/2011	16/07/2015	16/07/2015	
Initials	JG	TT	PJ	

Material			
No.	Plan	Article	Quantity
1	X004688	1734104	2
2	M10 DIN125A	1001558	2
3	Ø32	1001610	2
4	---	1904114	2
5	---	1910944	1
6	---	1910945	1
7	M8x16	1991858	2
8	M8 DIN934	1001671	2
9	M10x1 DIN934	1901345	2
10	M8x20 DIN7991/10.9	1904812	8
11	Block gray	2990480	2
	Line gray	1908436 (L)	1
	Line gray	1908437 (R)	1

E000566 / F

MOTOR 6KM/h SQUOD / SQOUD 2009 / SQOUD SU				Plan No.
				E000566 / F
Scale: 1:5	Drafted	Verified	Confirmed	
Date	11/02/2011	16/07/2015	16/07/2015	
Initials	JG	TT	PJ	