


Use WEICONLOCK AN 302-50

Use WEICONLOCK AN 302-43

Detail B  
Scale: 3:5

Drawing shows only right side

BACKREST ADJUSTABLE ANGLE TRIGO				Plan No.
				EP00356/ E
Scale: 3:10	Drafted	Verified	Confirmed	
Date	29/10/2018	03/06/2020	03/06/2020	
Initials	PM	MK	TK	

Material			
No.	Plan	Article	Quantity
* 1	YP06081(35)	2793265	1
	YP06082(37)	2793365	1
	YP06083(39)	2793465	1
	YP06084(41)	2793565	1
	YP06085(43)	2793665	1
	YP06086(45)	2793765	1
	YP06087(47)	2793865	1
	YP06088(49)	2793965	1
	YP06089(51)	2794065	1
2	XP01632	GP09399	1
3	XP02061	GP14700	1
4	XP01862	GP13200	1
5	M4x12	2994126	2
6	---	---	---
7	XP01861	GP12900	1
8	M5x12 DIN 7991	2994076	2
9	CAP $\phi$ 22	1990648	1
** 10	Fold. push handle	EP00246	1
** 11	M6x10	2994127	4
12	YP08333(35)	2713865	1
	YP08334(37)	2713965	1
	YP08335(39)	2714065	1
	YP08336(41)	2714165	1
	YP08337(43)	2714265	1
	YP08338(45)	2714365	1
	YP08339(47)	2714465	1
	YP08340(49)	2714565	1
	YP08341(51)	2714665	1
13	Hand grip GR-KS	2993809	1
	GR-ERGO-S	2993810	1

\* Gdy uzywamy z raczka skladana(10,11) wybrac 2 rozmiary krotsza rure.

\* When using a folded handle (10,11), select 2 sizes shorter tube.

\*\* Option.

\*\* Opcja.

EP00356 / E

BACKREST ADJUSTABLE ANGLE TRIGO				Plan No.
Page 2/2				EP00356/ E
Scale: 3:10	Drafted	Verified	Confirmed	
Date	29/10/2018	03/06/2020	03/06/2020	
Initials	PM	MK	TK	