



Use WEICONLOCK AN 302-50.

Use WEICONLOCK AN 302-43.

Detail B
Scale: 3:5

Drawing shows only right side

EP00374 / C

BACKREST ADJUSTABLE ANGLE AND ADJUSTABLE HIGHT WITH FOLDING HANDLE TRIGO				Plan No. EP00374 / C
Scale: 3:10	Drafted	Verified	Confirmed	
Date	13/02/2019	20/02/2020	20/02/2020	
Initials	AS	FD	TK	

Material			
No.	Plan	Article	Quantity
1	YP07083(35-45)	2745665	1
	YP07084(40-50)	2745765	1
2	YP07062(32-37)	2744565	1
	YP07060(35-40)	2744665	1
	YP07061 (40-45)	2744765	1
	YP07063(45-50)	2744865	1
3	XP02061	GP14700	1
4	XP01862	GP13200	1
5	M4x12	2994126	1
6	XP01632	GP09399	2
7	XP01861	GP12900	1
8	M5x12 DIN 7991	2994076	1
9	M6x30	1992721	2
10	M6	2994059	1
11	WASHER_M6	2994068	1
12	Fold. push handle	EP00246	1
13	Y004353	F000411	1
14	M6x10 DIN916	2994284	2
15	M6X10	2994127	2
*16	Cap D22	1990648	1

* When unused folding handle / Gdy nie uzywana raczka skladana.

EP00374 / C

BACKREST ADJUSTABLE ANGLE AND ADJUSTABLE HIGHT WITH FOLDING HANDLE TRIGOPage 2/2				Plan No. EP00374/ C
Scale: 3:10	Drafted	Verified	Confirmed	
Date	13/02/2019	20/02/2020	20/02/2020	
Initials	AS	FD	TK	